



**BHT**<sup>®</sup>  
Hygienetechnik

**Atherton**

# Load Carrier Systems

for the medical field

*Flexible systems for a wide range of wash goods*



# Series 3

## RACK FOR 6 INSTRUMENT TRAYS

**71080.06** ..... For 6 DIN-instrument trays, dimensions 485 x 250 mm, height 50 or 70 mm, 3 levels with 2 integrated rotating arms

■ **Following trays can be used in this rack:**  
 92090.01 ..... DIN-instrument tray 1/1 (485x250x50 mm)  
 92090.02 ..... DIN-instrument tray 1/2 (250x250x50 mm)  
 92100.01 ..... Tray for microsurgical instruments 1/1 (485x250x50 mm)  
 92100.02 ..... Tray for microsurgical instruments 1/2 (250x250x50 mm)

■ **Max. loading height per level:**  
 lowest level 70 mm, central level 100 mm, upper level up to max. 150 mm



## RACK FOR 8 INSTRUMENT TRAYS

**71080.08** ..... For 8 DIN-instrument trays, dimensions 485 x 250 mm, height 50 mm, 4 levels with 3 integrated rotating arms

■ **Following trays can be used in this rack:**  
 92090.01 ..... DIN-instrument tray 1/1 (485x250x50 mm)  
 92090.02 ..... DIN-instrument tray 1/2 (250x250x50 mm)  
 92100.01 ..... Tray for microsurgical instruments 1/1 (485x250x50 mm)  
 92100.02 ..... Tray for microsurgical instruments 1/2 (250x250x50 mm)

■ **Max. loading height per level:**  
 50 mm



## MIS-RACK

**71070.01** ..... For rigid endoscopes and instruments for minimal-invasive surgery; for 2 complete MIS-sets. With patented pressure chamber system.

**71070.02** ..... With half pressure chamber for urology and bronchoscopy



## SPIRAL RACK FOR ANAESTHETIC UTENSILS

- 7110.03 ..... For:
- 9 respiratory tubes up to 1,5 m length
  - 2 respiratory bags
  - 28 tubes and catheters
  - 4 respiratory masks
  - 2 secretion glasses
  - 1 corrugated tube
  - 1 circular item



## UNIVERSAL RACK

- 7120.00 ..... 2 levels, with rotating jet arm, for big wash goods, theatre shoes, kidney bowls, etc.

■ **Following baskets can be used in the universal rack:**

- 72120.03 ..... Basket for kidney bowls for upper and lower level for 12 kidney bowls



## BASKET FOR THEATRE SHOES

- 71125.00 ..... 2 level, for total 24 theatre shoes

## BASIC CART LOWER LEVEL

- 71010.02 ..... For mass goods, big wash goods, etc., or other wash goods in baskets



# Series 3

## BASIC CART UPPER LEVEL

71010.01 ..... For mass goods, big wash goods, etc.

■ **Following accessories are necessary:**  
Telescope rails for the use of upper level

51010.00.....



## DIRECT INJECTION RACK WITH 16 JETS

71040.16 ..... For suction bottles with max. Ø 112 mm  
Standard jet 6 x 100 mm



## RACK FOR BABY BOTTLES WITH 49 JETS

71050.49 (M3) ... Capacity: 1 basket for 49 baby bottles

■ **Following baskets can be used:**

Basket for 49 baby bottles,

92080.01 ..... content 250 ml max. OD

Basket for 49 baby bottles,

92080.03..... content 90 – 125 ml max. OD



## RACK FOR BABY BOTTLES WITH 42 JETS

71050.42..... Capacity: 2 baskets with 21 baby bottles each

■ **Following baskets can be used:**

Basket for 21 baby bottles,

92080.02..... content 250 ml max. OD

Basket for 21 baby bottles,

92080.04..... content 90 – 125 ml max. OD



# 3-4

Series 3

## MAGNETIC FIXTURES FOR RACKS

51008.00 . . . . . Necessary for each rack, if the option „Rack coding”  
#51007.00 is ordered

## Racks for flexible endoscopes and accessories

### RACKS FOR 2 FLEXIBLE ENDOSCOPES

20000.44 (E3) . . . With 2 independent pressure chambers, placed above  
each other, front loading

20000.45 (E4) . . . 1 basket for small items with lid and one blind plug for the  
pressure chamber included. Space for one secretion bottle.  
Docking facility for jet- resp. rinse channel adaptor

#### ■ Additional accessories:

Hose adapter for Olympus Exera 160 series

Hose adapter for Pentax 70K and 80K endoscopes

20000.07 . . . . . Basket for angioscopes

Reducing basket for smaller endoscopes, children's endoscopes  
or bronchoscopes with thinner supply tube and smaller exami-  
nation part; to be used in the pressure chamber.



# Series 3-4

Series 3

## ENDO-SONO-RACK

**20000.02 (E3) ...** For combined cleaning: on the left side, a big pressure chamber for an endo-sonoscope, on the right side a normal pressure chamber for a conventional endoscope.




## RACK FOR FLEXIBLE ACCESSORIES

**20000.10 (E3) ...** ■ For the reception of:  
 PE-pliers  
 ERCP-catheters  
**20000.42 (E4) ...** foreign body forceps  
 Papillotoms  
 Sclerosir needle  
 extractor baskets  
 Polypectomy wire  
 procto-rectoscopes  
 rinse bottles of all kinds



## BASKET WITH LID FOR SMALL ITEMS

**92100.40 ...** For use in the endoscope rack  
 dimensions 220 x 140 x 50 mm



## BLIND PLUG

**16391.00 ...** For use in the pressure chamber,  
 if only one endoscope is processed

# 4-5

Series 4-5

## RACK FOR 8 INSTRUMENT TRAYS

- 11080.08 (M4) ..** For 8 DIN-instrument trays, dimensions 485 x 250 mm, up to 100 mm height, 4 levels with 3 integrated rotating arms
- 11080.08S ..... (M5)**
- 92090.01 .....** Tray for microsurgical instruments 1/1 (485x250x50 mm)  
Tray for microsurgical instruments 1/2 (250x250x50 mm)

- **Following trays can be used in this rack:**  
DIN-instrument tray 1/1 (485x250x50 mm)  
DIN-instrument tray 1/2 (250x250x50 mm)  
Tray for microsurgical instruments 1/1 (485x250x50 mm)  
Tray for microsurgical instruments 1/2 (250x250x50 mm)
- **Max. loading height per level:**  
M4: max. 100 mm, upper level max. 80 mm, M5: max. 100 mm



## RACK FOR 10 INSTRUMENT TRAYS

- 11080.10 (M4) ..** For 10 DIN-instrument trays, dimensions 485 x 250 mm, 5 levels with 4 integrated rotating arms
- 11080.10S (M5) ..**
- 92090.01 .....** Tray for microsurgical instruments 1/1 (485x250x50 mm)  
Tray for microsurgical instruments 1/2 (250x250x50 mm)

- **Following trays can be used in this rack:**  
DIN-instrument tray 1/1 (485x250x50 mm)  
DIN-instrument tray 1/2 (250x250x50 mm)  
Tray for microsurgical instruments 1/1 (485x250x50 mm)  
Tray for microsurgical instruments 1/2 (250x250x50 mm)
- **Max. loading height per level:**  
M4: max. 65 mm; M5: max. 65 mm



## RACK FOR 12 INSTRUMENT TRAYS

- 11080.08T (M5-ISO) (T5-ISO)** For 12 DIN-instrument trays, dimensions 485 x 250 mm, height up to 100 mm, or 8 ISO-instrument trays, dimensions 485 x 344 x 50 mm, 4 levels with 6 integrated rotating arms
- 92090.01 .....** DIN-instrument tray 1/1 (485x250x50 mm)  
**92090.02 .....** DIN-instrument tray 1/2 (250x250x50 mm)  
**92100.01 .....** Tray for microsurgical instruments 1/1 (485x250x50 mm)  
**92100.02 .....** Tray for microsurgical instruments 1/2 (250x250x50 mm)  
**92090.01l .....** ISO-instrument tray 1/1 (485x344x50 mm)  
**92100.01l .....** ISO-tray for microsurgical instruments 1/1 (485x344x50 mm)

- **Max. loading height per level:** max. 100 mm



# 4-5

## Series 4-5

### RACK FOR 15 INSTRUMENT TRAYS

**11080.10T (M5-ISO)**  
**(T5-ISO)** For 15 DIN-instrument trays,  
dimensions 485 x 250 mm, height 50 mm,  
or 10 ISO-instrument trays,  
dimensions 485 x 344 x 50 mm  
5 levels with 8 integrated rotating arms

92090.01  
92090.02.....  
92100.01.....  
92100.02.....  
92090.01l.....  
92100.01l.....

■ **Following trays can be used in this rack:**  
DIN-instrument tray 1/1 (485x250x50 mm)  
DIN-instrument tray 1/2 (250x250x50 mm)  
Tray for microsurgical instruments 1/1 (485x250x50 mm)  
Tray for microsurgical instruments 1/2 (250x250x50 mm)  
ISO-instrument tray 1/1 (485x344x50 mm)  
ISO-tray for microsurgical instruments 1/1  
(485x344x50 mm)

■ **Max. loading height per level:**  
max. 50 mm



### MIS-RACK

**11070.01 (M4)**  
**11070.01S (M5)**  
**11070.01T (M5-ISO)**

For rigid endoscopes and instruments for  
minimal-invasive surgery; for 3 complete MIS-sets  
(ISO for 3/4 complete MIS-sets).  
With patented pressure chamber system.



### ROTO-SYSTEM

**11100.00 (M4)...**  
**11100.00S (M5)...**  
**11100.00l (M5-ISO)**

Rotating hose rinsing appliance for 11 AN-sets  
(M5 + ISO: 13 sets), = 22 (26) corrugated tubes,  
14 (17) suction tubes and further accessories.  
The rotation of the drum opposite to the flow direction  
of the water is the only possible solution to clean and  
disinfect the interior of the corrugated respiratory tubes.  
The even distribution of the residual water during the  
drying facilitates a fast and efficient drying.

- Hose magazine and covering net included -



# 4-5

Series 4-5

## SPIRAL RACK FOR ANAESTHETIC UTENSILS

1110.03 (M4)...	■ Capacity:	(ISO)
1110.03S (M5)...	12	(M4), 16 (M5)
1110.03T (M5-ISO)	respiratory tubes up to 1,5 m length	(16)
(T5-ISO)	6 respiratory bags	(9)
	16 tubes and catheters	(22)
	5 respiratory masks	(7)
	3 secretion glasses	(5)
	2 corrugated tubes	(3)
	1 circular item	(2)
	4 Y-pieces	(6)



## RACK FOR ANAESTHETIC ACCESSORIES

1110.04 (M4)...	For all systems		
1110.04S (M5)...	■ Capacity:	(ISO)	
1110.04T (M5-ISO)	4	respiratory tubes	(6)
(T5-ISO)	4 circular items		(6)
	4 respiratory chalk bags		(6)
	6 secretion bottles		(9)
	10 respiratory bags		(15)
	12 endotracheal tubes		(18)
	24 Guedel tubes		(36)
	with upper grid for small items etc.		



## DIRECT INJECTION RACK WITH 16 JETS

11040.16 (M4)...	For suction bottles with max. Ø 115 mm
	Standard jet 8 x 150 mm

## DIRECT INJECTION RACK WITH 24 JETS

11040.24S (M5)...	For suction bottles with max. Ø 100 mm
11040.24T (M5-ISO)	Standard jet 6 x 250 mm
(T5-ISO)	



# 4-5

## Series 4-5

### BASIC RACK

- 11010.01 (M4) . . . For mass goods, big wash goods, etc.  
11010.01S (M5) . . . to be used in the lower level  
11010.01T (M5-ISO)  
(T5-ISO)

■ For the use in the upper level (only M4/M5), following accessories are necessary:

- 11010.02 (M4) . . . Rotating jet arm and rails for the upper level,  
14390.00 (M5) . . . for the cleaning at 2 levels



### UNIVERSAL RACK

- 1120.00 (M4) . . . 2 levels, with rotating jet arm,  
1120.00S (M5) . . . for big wash goods, etc.  
1120.00T (M5-ISO)  
(T5-ISO)

■ Following baskets can be used in the universal rack:

- 12120.03 (M4) . . . Basket for kidney bowls for upper and lower level  
12120.03S (M5) . . . for 30 kidney bowls (250x155x44 mm)  
12120.03T (M5-ISO) **ISO:** for 39 kidney bowls  
(T5-ISO)
- 12120.04 (M4) . . . Basket for kidney bowls for upper and lower level  
12120.04S (M5) . . . for 12 kidney bowls or care sets  
12120.04T (M5-ISO) **M5:** 18 kidney bowls, **ISO:** 24 kidney bowls  
(T5-ISO)
- 12120.05 (M4) . . . Basket for bowls for lower level  
12120.05S (M5) . . . for 4 bowls (content 6 l)  
12120.05T (M5-ISO) **ISO:** for 6 bowls  
(T5-ISO)



### RACK FOR THEATRE SHOES

- 1121.00 (M4) . . . for 24 pairs of theatre shoes  
1121.00S (M5) . . . for 30 pairs of theatre shoes  
1121.00T (M5-ISO)  
(T5-ISO) for 40 pairs of theatres shoes

# 4-5

Series 4-5

## RACK FOR AESCULAP- / MARTIN/WAGNER CONTAINERS

- 11060.01 (M4) ... For 4 containers with lid or 6 containers with 2 lids  
dimensions 600 x 300 x 150 mm
- 11060.01S (M5) ...
- 11060.01T (M5-ISO) ISO: for 6 DIN-containers and 6 lids or  
4 ISO-containers with lid,  
dimensions 600 x 400 x 200 mm
- (T5-ISO)



## RACK FOR WAGNER CONTAINERS

- 11060.02 (M4) ... For 4 containers with lid,  
dimensions 600 x 300 x 150 mm



## RACK FOR 1 STERILE UNIT

- 11060.03 (M4) ... For 2 containers with lid, dimensions 600 x 300 x 300 mm
- 11060.03S (M5) ... For 3 containers with lid, dimensions 600 x 300 x 300 mm
- 11060.03T (M5-ISO) 4 ISO-containers with lid, dimensions 600 x 400 x 200 mm  
(T5-ISO) or 6 containers with lid, dimensions 600 x 300 x 150 mm



## COMBI RACK

- 11060.05 (M4) ... For 4 DIN-instrument trays, 485 x 250 x 50 mm,  
and 2 containers, DIN-dimensions 600 x 300 x 150 mm,  
with lid (only for containers type Aesculap / Martin)
- 11060.05S (M5) ...
- 11060.05T (M5-ISO) ISO: for 6 instrument trays and 2 containers  
(T5-ISO)



# 4-5

## Series 4-5

### RACK FOR BABY BOTTLES WITH 63 JETS

**11050.63 (M4)**... Capacity: 3 baskets with 21 baby bottles each

■ **Following baskets can be used:**

- 92080.02 ..... Basket for 21 baby bottles, content 250 ml
- 92080.04 ..... Basket for 21 baby bottles, content 90 – 125 ml
- 92080.05 ..... Basket for teats (2pcs.)



### RACK FOR BABY BOTTLES WITH 70/84 JETS

**11050.70 (M4)**... 1 basket with 49 + 1 basket  
**11050.70S (M5)**... with 21 baby bottles

■ **Capacity:**

**11050.70T (M5-ISO)** 4 baskets with each 21 baby bottles  
**11050.84T (T5)**... Add. upper level with jet arm for teats

■ **Following baskets can be used:**

- Basket for 49 baby bottles, content 250 ml
- 92080.01 ..... Basket for 49 baby bottles, content 90 – 125 ml
- 92080.03 ..... Basket for 21 baby bottles, content 250 ml
- 92080.02 ..... Basket for 21 baby bottles, content 90 – 125 ml
- 92080.04 ..... Basket for teats (3 pcs.)
- 92080.05 .....



### SPECIAL RACK FOR BABY BOTTLES

**11050.00 (M4)**... 2-level rack with 63 jets at the lower level;  
**11050.00S(M5)**... upper level with rotating jet arm.  
For 63 baby bottles with teats

■ **Following baskets can be used:**

- 92080.02 or ..... (3x) Basket for 21 baby bottles, content 250 ml
- 92080.04 ..... (3x) Basket for 21 baby bottles, content 90 – 125 ml
- 92080.05 and ... (2x) Basket for teats

### MAGNETIC FIXTURES FOR RACKS

**52008.00(M4)** ... Necessary for each rack, if the option „Rack coding“  
**52009.00(M5)** ... # 52007.00/52007.10 is ordered  
(M5-ISO)

**P: 1300 580 870**

**E: salesandservice@atherton.net**

**W: www.atherton.net**

**TDS 140 Revision 0 13/12/2010**



Certified by DQS according to  
DIN EN ISO 9001: 2000  
DIN EN 13485: 2003  
Reg. No. 019906-01

BHT Hygienetechnik GmbH  
Messerschmittstraße 11  
D - 86368 Gersthofen  
Tel. +49-(0)8 21/2 78 93-0  
Fax +49-(0)8 21/78 40 99  
E-Mail: hygiene@bht.de  
Internet: www.bht.de

# Atherton